

[AFYS]3G

Laser marks directly on the surface of tubes

Laser any shape, text, logo or sample ID on a tube



INFORMATION MARKING SYSTEM

The AFYS3G Manual Information Marking System Lambda8 can mark tubes from 0.30ml to 15.00ml. Multiple adapters are available to accommodate the wide variety of tubes with which the Lambda8 is compatible. With this device, no labels, stickers or ink is used to mark tubes. The laser engraves markings with any shape or text directly onto the surface of each tube. The markings are resistant against chemicals, mechanical abrasion and temperatures between +100°C and -196°C, ensuring sample traceability.

ARTICLE #	DESCRIPTION	DIMENSIONS (MM)
A20101	AFYS3G Manual Information Marking System Lambda8	482 x 416 x 778
A20102	AFYS3G Manual Information Marking System Lambda8 + Integrated PC	482 x 416 x 778*
A20501	AFYS3G Air Filtration System	420 x 425 x 400

*PC not included