

# CappAero384

MULTICHANNEL PIPETTES

## CappAero384™ Multichannel Pipettes for 384-well plates

Most applications involving 384-well plates are performed using automated systems. Capp provides the world's only alternative for the quick and high precision pipetting of 16-, 48-, or 64-wells at once without the use of expensive, automated equipment.

## Ideal for 384 well liquid handling

The CappAero384™ 48 & 64 channel pipette allows for pipetting large samples within a short period of time while reducing the workload. The CappAero384™ is generally applicable for all transfers of liquids from reservoir and 384 deep-well plates into 384 well test plates, and it offers very obvious advantages in PCR, micro array, and MALDI spectrometry techniques. CappAero384™ Multichannel Pipettes feature lightweight construction and produce a consistent and reliable liquid transfer result in one ergonomic solution.



## ExpellPlus384™ low retention

## High quality, low retention tips

Capp offers a range of ExpellPlus384™ tips which are specifically designed to perfectly fit the Capp range of 384 well pipetting. Due to special polymers and an extremely smooth inside surface, the ExpellPlus384™ tips from Capp virtually eliminate sample retention, and reduce the binding of DNA, RNA, and Proteins to a minimum.

The range also includes a filter tip version to effectively prevent cross contamination due to aerosol formation (such as DNA), and dangers from toxic, infectious, radioactive or other hazardous material's carry over.



# Ordering Information

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
C02-16	Capp multi pipettes, 16-channel, 0.2-2µl	4.00	2.00	●
C10-16	Capp multi pipettes, 16-channel, 0.5-10µl	2.00	1.00	○
C50-16	Capp multi pipettes, 16-channel, 5-50µl	1.40	0.60	●
C10-48	Capp multi pipettes, 48-channel, 0.5-10µl	4.00	2.00	○
C10-64	Capp multi pipettes, 64-channel, 0.5-10µl	4.00	2.00	○

## ExpellPlus384™ low retention tips

Cat. No.	Description	Qty/pack	Packs/case
5030005C	ExpellPlus 0.2-50µl, racked, sterile	10 x 384	5
5030006C	ExpellPlus 0.2-50µl, racked, sterile, filtered	10 x 384	5

## Expell384 PCR plates

Cat. No.	Description
5100513C	Expell 384-well PCR plate, natural, 320 pcs.



More information on page 38.